

SDK version: 2017.6.6.1

\* Call the NIREz20\_Initial () before open the SDK function

----- SDK function -----

1. NIREz20 initialization:

Int NIREz20\_Initial (void);

Return value: 1 is normal

Parameters: None

2. Connection:

Void NIREz20\_Connect (void);

Return value: none

Parameters: None

3. NIREz20 Link Status:

Int NIREz20\_Link\_Status (void);

Return value:

0 -> is not connected

1 -> is connected

-99 -> incomplete initialization

Parameters: None

4. Write scan parameters:

int NIREz20\_Config ( int i\_lamp\_switch , int i\_start\_wavelength , int  
i\_end\_wavelength , double d\_width , int i\_gain ,int i\_points ,  
int i\_average,int i\_exposure,int i\_scan\_mode );

Return value:

0 -> OK

-1 -> has not yet read the software license

-2 -> USB is not connected to NIREz

-3 -> LAMP parameter setting error

-4 -> Start wavelength parameter setting error

-5 -> End Wavelength parameter setting error

-6 -> Wavelength Width parameter setting error

-7 -> Gain parameter setting error

-8 -> Points parameter setting error

-9 -> The mean parameter setting is incorrect

-10 -> Exposure time parameter setting error

-11 -> The scan mode is set incorrectly

-99 -> incomplete initialization

Parameters:

Int i\_lamp\_switch: 0 is LAMP OFF; 1 is LAMP ON

Int i\_start\_wavelength: Start wavelength 900 ~ 1700 nm

Int i\_end\_wavelength: End wavelength 900 ~ 1700 nm

Double d\_width: Wavelength Width (2.34 to 60.84 nm) The accumulated value is 1.17

Int i\_gain: gain 1, 2, 4, 8, 16, 32,

Int i\_points: the number of measured points

Int i\_average: Set the average times after scanning (1 to 99)

Int i\_exposure: Set the exposure time value (0 ~ 6)

0 -> 0.635ms

1 -> 1.27ms

2 -> 2.45ms

3 -> 5.08ms

4 -> 15.24ms

5 -> 30.48ms

6 -> 60.96ms

Int i\_scan\_mode: Set the scan mode

0 -> Column mode

1 -> Hadamard mode

5. Execute a scan:

Int Nirez20\_Scan (int \* i\_total, int \* i\_gain, double \* d\_wavelength, int \* i\_intensity);

Front Return Value: Returns the scan status

0 -> OK

-1 -> has not yet read the software license

-2 -> USB is not connected to NIREz

-3 -> Software authorization code error

-4 -> The scan parameters are not configured first

-5 -> scan is in progress

-6 -> scan results are abnormal

-7 -> scan results overflow

-99 -> incomplete initialization

Parameter Return Value:

Int \* i\_total: The number of data

Int \* i\_gain: gain

Double \* d\_wavelength: Wavelength array

Int \* i\_intensity: Intensity array

Note: scan results, without any treatment

The license key will be checked before scanning

6. Specify the License Key file path:

Void NIRez20\_Set\_Key\_Path (char \* c\_path);

Return value: none

Parameters:

Char \* c\_path: set the path where the License Key file is placed (Example -> D: \ NIR\_KEY)

Currently, the license key file has two paths:

1. C: \ Isuzu Optics \ NIRez \ NIRez\_License\_Key \ SDK
2. and NIRez20\_SDK.dll under the same path

If the above two paths can not find the License Key file, use this function to specify the path

\* If the License Key file is not found, it will return -3 when NIRez20\_Scan () is executed

7. Read serial number:

Void NIRez20\_Get\_Sn (char \* c\_file, char \* c\_nirez);

Return value: none

Parameter Return Value:

Char \* c\_file: File not exists! If no License Key file is found; otherwise, return the serial number in the License Key file

Char \* c\_nirez: Returns the hardware number of NIRez